

ANKMY1 Antibody



PACO42030

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

ANKMY1

Uniprot

Q9P2S6

Synonyms:

Ankyrin repeat and MYND domain-containing protein 1 (Testis-specific ankyrin-like protein 1) (Zinc finger MYND domain-containing protein 13), ANKMY1, TSAL1 ZMYND13


Immunogen:

Recombinant Human Ankyrin repeat and MYND domain-containing protein 1 protein (1-220AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO42030 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO42030 at dilution of 1:100.