

REEP2 Antibody, Biotin conjugated



PACO42093

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for endoplasmic reticulum (ER) network formation, shaping and remodeling. May enhance the cell surface expression of odorant receptors.

Gene ID:

REEP2

Uniprot

Q9BRK0

Synonyms:

Receptor expression-enhancing protein 2, REEP2, C5orf19 SGC32445

Immunogen:

Recombinant Human Receptor expression-enhancing protein 2 protein (78-252AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A