

MECR Antibody, HRP conjugated



PACO42107

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Oxidoreductase with a preference for short and medium chain substrates, including trans-2-hexenoyl-CoA (C6), trans-2-decenoyl-CoA (C10), and trans-2-hexadecenoyl-CoA (C16). May play a role in mitochondrial fatty acid, synthesis.

Gene ID:

MECR

Uniprot

Q9BV79

Synonyms:

Enoyl-[acyl-carrier-protein] reductase, mitochondrial (EC 1.3.1.104) (2-enoyl thioester reductase) (Nuclear receptor-binding factor 1) (HsNrnf-1) (NRBF-1), MECR, NBRF1

Immunogen:

Recombinant Human Trans-2-enoyl-CoA reductase, mitochondrial protein (211-373AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A