

# NUP58 Antibody, HRP conjugated



PACO42115

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Component of the nuclear pore complex, a complex required for the trafficking across the nuclear membrane.

**Gene ID:**

NUP58

**Uniprot**

Q9BVL2

**Synonyms:**

Nucleoporin p58/p45 (58 kDa nucleoporin) (Nucleoporin-like protein 1), NUP58, KIAA0410 NUPL1

**Immunogen:**

Recombinant Human Nucleoporin p58/p45 protein (211-480AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A