

FGFBP2 Antibody



PACO42138

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

Gene ID:

FGFBP2

Uniprot

Q9BYJ0

Synonyms:

Fibroblast growth factor-binding protein 2 (FGF-BP2) (FGF-binding protein 2) (FGFBP-2) (37 kDa killer-specific secretory protein) (Ksp37) (HBp17-related protein) (HBp17-RP), FGFBP2, KSP37


Immunogen:

Recombinant Human Fibroblast growth factor-binding protein 2 protein (20-223AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: FGFBP2 antibody at 2 μ g/ml + Human serum. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 25 kDa. Observed band size: 25 kDa.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO42138 at dilution of 1:100.