

ISCU Antibody, Biotin conjugated



PACO42157

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the assembly or repair of the [Fe-S] clusters present in iron-sulfur proteins. Binds iron.

Gene ID:

ISCU

Uniprot

Q9H1K1

Synonyms:

Iron-sulfur cluster assembly enzyme ISCU, mitochondrial (NifU-like N-terminal domain-containing protein) (NifU-like protein), ISCU, NIFUN

Immunogen:

Recombinant Human Iron-sulfur cluster assembly enzyme ISCU, mitochondrial protein (35-167AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A