

PORCN Antibody, Biotin conjugated



PACO42161

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Protein-serine O-palmitoleyltransferase that acts as a key regulator of the Wnt signaling pathway by mediating the attachment of palmitoleate, a 16-carbon monounsaturated fatty acid, (C16:1), to Wnt proteins. Serine palmitoleylation of WNT proteins is required for efficient binding to frizzled receptors.

Gene ID:

PORCN

Uniprot

Q9H237

Synonyms:

Protein-serine O-palmitoleyltransferase porcupine (EC 2.3.1. -) (Protein MG61), PORCN, MG61 PORC PPN

Immunogen:

Recombinant Human Protein-serine O-palmitoleyltransferase porcupine protein (269-332AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A