

# ARRDC1-AS1 Antibody, Biotin conjugated



PACO42165

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

ARRDC1-AS1

**Uniprot**

Q9H2J1

**Synonyms:**

Uncharacterized protein ARRDC1-AS1 (ARRDC1 antisense RNA 1) (ARRDC1 antisense gene protein 1), ARRDC1-AS1, C9orf37

**Immunogen:**

Recombinant Human Uncharacterized protein ARRDC1-AS1 protein (1-176AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A