

# ZCCHC4 Antibody, Biotin conjugated



PACO42173

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be a methyltransferase.

**Gene ID:**

ZCCHC4

**Uniprot**

Q9H5U6

**Synonyms:**

Zinc finger CCHC domain-containing protein 4, ZCCHC4

**Immunogen:**

Recombinant Human Zinc finger CCHC domain-containing protein 4 protein (25-229AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A