

RPP21 Antibody, HRP conjugated



PACO42179

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5'-ends.

Gene ID:

RPP21

Uniprot

Q9H633

Synonyms:

Ribonuclease P protein subunit p21 (RNaseP protein p21) (EC 3.1.26.5) (Ribonuclease P/MRP 21 kDa subunit) (Ribonucleoprotein V), RPP21, C6orf135 CAT60

Immunogen:

Recombinant Human Ribonuclease P protein subunit p21 protein (2-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A