

# RPP21 Antibody, FITC conjugated



PACO42180

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Component of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5'-ends.

**Gene ID:**

RPP21

**Uniprot**

Q9H633

**Synonyms:**

Ribonuclease P protein subunit p21 (RNaseP protein p21) (EC 3.1.26.5) (Ribonuclease P/MRP 21 kDa subunit) (Ribonucleoprotein V), RPP21, C6orf135 CAT60

**Immunogen:**

Recombinant Human Ribonuclease P protein subunit p21 protein (2-154AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A