

ADAP2 Antibody, HRP conjugated



PACO42219

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

GTPase-activating protein for the ADP ribosylation factor family (Potential). Binds phosphatidylinositol 3,4,5-trisphosphate (PtdInsP3) and inositol 1,3,4,5-tetrakisphosphate (InsP4). Possesses a stoichiometry of two binding sites for InsP4 with identical affinity.

Gene ID:

ADAP2

Uniprot

Q9NPF8

Synonyms:

Arf-GAP with dual PH domain-containing protein 2 (Centaurin-alpha-2) (Cnt-a2), ADAP2, CENTA2

Immunogen:

Recombinant Human Arf-GAP with dual PH domain-containing protein 2 protein (110-381AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A