

MANBAL Antibody, Biotin conjugated



PACO42225

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	MANBAL
Source:	Uniprot
Rabbit	Q9NQG1
Isotype:	Synonyms:
IgG	Protein MANBAL, MANBAL
Applications:	Immunogen:
ELISA	Recombinant Human Protein MANBAL protein (45-85AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A