

EXOSC3 Antibody, HRP conjugated



PACO42227

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EXOSC3

Uniprot

Q9NQT5

Synonyms:

Exosome complex component RRP40 (Exosome component 3) (Ribosomal RNA-processing protein 40) (p10), EXOSC3, RRP40

Immunogen:

Recombinant Human Exosome complex component RRP40 protein (2-275AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A