

# NLE1 Antibody, HRP conjugated



PACO42267

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Plays a role in regulating Notch activity. Plays a role in regulating the expression of CDKN1A and several members of the Wnt pathway, probably via its effects on Notch activity. Required during embryogenesis for inner mass cell survival.

**Gene ID:**

NLE1

**Uniprot**

Q9NVX2

**Synonyms:**

Notchless protein homolog 1, NLE1

**Immunogen:**

Recombinant Human Notchless protein homolog 1 protein (2-485AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A