

NLE1 Antibody, FITC conjugated



PACO42268

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in regulating Notch activity. Plays a role in regulating the expression of CDKN1A and several members of the Wnt pathway, probably via its effects on Notch activity. Required during embryogenesis for inner mass cell survival.

Gene ID:

NLE1

Uniprot

Q9NVX2

Synonyms:

Notchless protein homolog 1, NLE1

Immunogen:

Recombinant Human Notchless protein homolog 1 protein (2-485AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A