

ACTR10 Antibody



PACO42290

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Gene ID:

ACTR10

Uniprot

Q9NZ32

Synonyms:

Actin-related protein 10 (Actin-related protein 11) (hARP11), ACTR10, ACTR11 ARP11


Immunogen:

Recombinant Human Actin-related protein 10 protein (165-417AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: PC-3 whole cell lysate, HEK293 whole cell lysate. All lanes: ACTR10 antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 47 kDa. Observed band size: 47 kDa.