

ACTR10 Antibody, Biotin conjugated



PACO42293

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ACTR10

Uniprot

Q9NZ32

Synonyms:

Actin-related protein 10 (Actin-related protein 11) (hARP11), ACTR10, ACTR11 ARP11

Immunogen:

Recombinant Human Actin-related protein 10 protein (165-417AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A