

JPT1 Antibody, FITC conjugated



PACO42332

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

JPT1

Uniprot

Q9UK76

Synonyms:

Jupiter microtubule associated homolog 1 (Androgen-regulated protein 2) (Hematological and neurological expressed 1 protein) [Cleaved into: Jupiter microtubule associated homolog 1, N-terminally processed], JPT1, ARM2 HN1

Immunogen:

Recombinant Human Jupiter microtubule associated homolog 1 protein (2-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A