

BAG5 Antibody, HRP conjugated



PACO42339

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibits both auto-ubiquitination of PARK2 and ubiquitination of target proteins by PARK2 By similarity. May function as a nucleotide exchange factor for HSP/HSP70, promoting ADP release, and activating Hsp70-mediated refolding.

Gene ID:

BAG5

Uniprot

Q9UL15

Synonyms:

BAG family molecular chaperone regulator 5 (BAG-5) (Bcl-2-associated athanogene 5), BAG5, KIAA0873

Immunogen:

Recombinant Human BAG family molecular chaperone regulator 5 protein (1-447AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A