

STAP1 Antibody



PACO42342

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

In BCR signaling, appears to function as a docking protein acting downstream of TEC and participates in a positive feedback loop by increasing the activity of TEC.

Gene ID:

STAP1

Uniprot

Q9ULZ2

Synonyms:

Signal-transducing adaptor protein 1 (STAP-1) (BCR downstream-signaling protein 1) (Docking protein BRDG1) (Stem cell adaptor protein 1), STAP1, BRDG1

Immunogen:

Recombinant Human Signal-transducing adaptor protein 1 protein (1-295AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A