

SNX6 Antibody, FITC conjugated



PACO42348

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SNX6

Uniprot

Q9UNH7

Synonyms:

Sorting nexin-6 (TRAF4-associated factor 2) [Cleaved into: Sorting nexin-6, N-terminally processed], SNX6

Immunogen:

Recombinant Human Sorting nexin-6 protein (1-290AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A