

# DSE Antibody, HRP conjugated



PACO42359

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Converts D-glucuronic acid, to L-iduronic acid, (IdoUA) residues.

**Gene ID:**

DSE

**Uniprot**

Q9UL01

**Synonyms:**

Dermatan-sulfate epimerase (DS epimerase) (EC 5.1.3.19) (Chondroitin-glucuronate 5-epimerase) (Squamous cell carcinoma antigen recognized by T-cells 2) (SART-2), DSE, SART2

**Immunogen:**

Recombinant Human Dermatan-sulfate epimerase protein (23-270AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A