

# ZNRF3 Antibody, Biotin conjugated



PACO42365

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

E3 ubiquitin-protein ligase that acts as a negative regulator of the Wnt signaling pathway by mediating the ubiquitination and subsequent degradation of Wnt receptor complex components Frizzled and LRP6. Acts on both canonical and non-canonical Wnt signaling pathway. Acts as a tumor suppressor in the intestinal stem cell zone by inhibiting the Wnt signaling pathway, thereby restricting the size of the intestinal stem cell zone.

**Gene ID:**

ZNRF3

**Uniprot**

Q9ULT6

**Synonyms:**

E3 ubiquitin-protein ligase ZNRF3 (EC 2.3.2.27) (RING finger protein 203) (RING-type E3 ubiquitin transferase ZNRF3) (Zinc/RING finger protein 3), ZNRF3, KIAA1133 RNF203

**Immunogen:**

Recombinant Human E3 ubiquitin-protein ligase ZNRF3 protein (56-219AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A