

ANGPTL2 Antibody, FITC conjugated



PACO42376

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Induces sprouting in endothelial cells through an autocrine and paracrine action.

Gene ID:

ANGPTL2

Uniprot

Q9UKU9

Synonyms:

Angiotensin-converting enzyme 2 (Angiotensin-converting enzyme 2), ANGPTL2, ARP2

Immunogen:

Recombinant Human Angiotensin-converting enzyme 2 protein (23-493AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A