

PADI1 Antibody, HRP conjugated



PACO42383

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the deimination of arginine residues of proteins.

Gene ID:

PADI1

Uniprot

Q9ULC6

Synonyms:

Protein-arginine deiminase type-1 (EC 3.5.3.15) (Peptidylarginine deiminase I) (Protein-arginine deiminase type I), PADI1, PAD1 PDI1

Immunogen:

Recombinant Human Protein-arginine deiminase type-1 protein (1-663AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A