

PLA2G2D Antibody, Biotin conjugated



PACO42401

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

PA2 catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides. L-alpha-1-palmitoyl-2-linoleoyl phosphatidylethanolamine is more efficiently hydrolyzed than the other phospholipids examined.

Gene ID:

PLA2G2D

Uniprot

Q9UNK4

Synonyms:

Group IID secretory phospholipase A2 (GIID sPLA2) (sPLA2-IIID) (EC 3.1.1.4) (PLA2IID) (Phosphatidylcholine 2-acylhydrolase 2D) (Secretory-type PLA, stroma-associated homolog), PLA2G2D, SPLASH

Immunogen:

Recombinant Human Group IID secretory phospholipase A2 protein (22-145AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A