

PSG11 Antibody, Biotin conjugated



PACO42409

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

female pregnancy.

Gene ID:

PSG11

Uniprot

Q9UQ72

Synonyms:

Pregnancy-specific beta-1-glycoprotein 11 (PS-beta-G-11) (PSBG-11) (Pregnancy-specific glycoprotein 11) (Pregnancy-specific beta-1-glycoprotein 13) (PS-beta-G-13) (PSBG-13) (Pregnancy-specific glycoprotein 13), PSG11, PSG13 PSG14

Immunogen:

Recombinant Human Pregnancy-specific beta-1-glycoprotein 11 protein (7-219AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A