

# VPREB3 Antibody, Biotin conjugated



PACO42417

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Associates with the Ig-mu chain to form a molecular complex that is expressed on the surface of pre-B-cells.

**Gene ID:**

VPREB3

**Uniprot**

Q9UK13

**Synonyms:**

Pre-B lymphocyte protein 3 (N27C7-2) (Protein VPreB3), VPREB3

**Immunogen:**

Recombinant Human Pre-B lymphocyte protein 3 protein (21-122AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A