

# FBXW2 Antibody, FITC conjugated



PACO42420

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

**Gene ID:**

FBXW2

**Uniprot**

Q9UKT8

**Synonyms:**

F-box/WD repeat-containing protein 2 (F-box and WD-40 domain-containing protein 2) (Protein MD6), FBXW2, FBW2 FWD2

**Immunogen:**

Recombinant Human F-box/WD repeat-containing protein 2 protein (72-340AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A