

FBXW2 Antibody, Biotin conjugated



PACO42421

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

Gene ID:

FBXW2

Uniprot

Q9UKT8

Synonyms:

F-box/WD repeat-containing protein 2 (F-box and WD-40 domain-containing protein 2) (Protein MD6), FBXW2, FBW2 FWD2

Immunogen:

Recombinant Human F-box/WD repeat-containing protein 2 protein (72-340AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A