

PACO42438

---

## Product Information

**Size:**

50ug

**Reactivity:**

Schizosaccharomyces pombe

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Phosphorylates both NADH and NAD(+), with a preference for NADH. Anti-oxidant factor and key source of the cellular reductant NADPH.

**Gene ID:**

pos5

**Uniprot**

Q9UT98

**Synonyms:**

NADH kinase pos5, mitochondrial (EC 2.7.1.86), pos5

**Immunogen:**

Recombinant Schizosaccharomyces pombe NADH kinase pos5, mitochondrial protein (43-386AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A