

SLC13A4 Antibody, Biotin conjugated



PACO42445

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Sodium/sulfate cotransporter that mediates sulfate reabsorption in the high endothelial venules (HEV).

Gene ID:

SLC13A4

Uniprot

Q9UKG4

Synonyms:

Solute carrier family 13 member 4 (Na(+)/sulfate cotransporter SUT-1) (NaS2), SLC13A4, SUT1

Immunogen:

Recombinant Human Solute carrier family 13 member 4 protein (134-273AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A