

ATXN7L1 Antibody, FITC conjugated



PACO42484

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	ATXN7L1
Source:	Uniprot
Rabbit	Q9ULK2
Isotype:	Synonyms:
IgG	Ataxin-7-like protein 1 (Ataxin-7-like protein 4), ATXN7L1, ATXN7L4 KIAA1218
Applications:	Immunogen:
ELISA	Recombinant Human Ataxin-7-like protein 1 protein (1-146AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A