

# ATXN7L1 Antibody, Biotin conjugated



PACO42485

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	ATXN7L1
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q9ULK2
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Ataxin-7-like protein 1 (Ataxin-7-like protein 4), ATXN7L1, ATXN7L4 KIAA1218
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human Ataxin-7-like protein 1 protein (1-146AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A