

# DIMT1 Antibody, FITC conjugated



PACO42492

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

pecifically dimethylates two adjacent adenosines in the loop of a conserved hairpin near the 3'-end of 18S rRNA in the 40S particle.

**Gene ID:**

DIMT1

**Uniprot**

Q9UNQ2

**Synonyms:**

Probable dimethyladenosine transferase (EC 2.1.1.183) (DIM1 dimethyladenosine transferase 1 homolog) (DIM1 dimethyladenosine transferase 1-like) (Probable 18S rRNA (adenine(1779)-N(6)/adenine(1780)-N(6))-dimethyltransferase) (Probable 18S rRNA dimethylase) (Probable S-adenosylmethionine-6-N', N'-adenosyl(rRNA) dimethyltransferase), DIMT1, DIMT1L

**Immunogen:**

Recombinant Human Probable dimethyladenosine transferase protein (1-313AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A