

DIMT1 Antibody, Biotin conjugated



PACO42493

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

pecifically dimethylates two adjacent adenosines in the loop of a conserved hairpin near the 3'-end of 18S rRNA in the 40S particle.

Gene ID:

DIMT1

Uniprot

Q9UNQ2

Synonyms:

Probable dimethyladenosine transferase (EC 2.1.1.183) (DIM1 dimethyladenosine transferase 1 homolog) (DIM1 dimethyladenosine transferase 1-like) (Probable 18S rRNA (adenine(1779)-N(6)/adenine(1780)-N(6))-dimethyltransferase) (Probable 18S rRNA dimethylase) (Probable S-adenosylmethionine-6-N', N'-adenosyl(rRNA) dimethyltransferase), DIMT1, DIMT1L

Immunogen:

Recombinant Human Probable dimethyladenosine transferase protein (1-313AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A