

# TRIM35 Antibody, Biotin conjugated



PACO42497

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role as a tumor suppressor. Implicated in the cell death mechanism.

**Gene ID:**

TRIM35

**Uniprot**

Q9UPQ4

**Synonyms:**

Tripartite motif-containing protein 35 (Hemopoietic lineage switch protein 5), TRIM35, HLS5 KIAA1098

**Immunogen:**

Recombinant Human Tripartite motif-containing protein 35 protein (1-206AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A