

PLCL2 Antibody



PACO42538

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May play an role in the regulation of Ins(1,4,5)P3 around the endoplasmic reticulum.

Gene ID:

PLCL2

Uniprot

Q9UPR0

Synonyms:

Inactive phospholipase C-like protein 2 (PLC-L(2)) (PLC-L2) (Phospholipase C-L2) (Phospholipase C-epsilon-2) (PLC-epsilon-2), PLCL2, KIAA1092 PLCE2


Immunogen:

Recombinant Human Inactive phospholipase C-like protein 2 protein (2-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO42538 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO42538 at dilution of 1:100.