

# PLCL2 Antibody, HRP conjugated



PACO42539

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May play an role in the regulation of Ins(1,4,5)P3 around the endoplasmic reticulum.

**Gene ID:**

PLCL2

**Uniprot**

Q9UPR0

**Synonyms:**

Inactive phospholipase C-like protein 2 (PLC-L(2)) (PLC-L2) (Phospholipase C-L2) (Phospholipase C-epsilon-2) (PLC-epsilon-2), PLCL2, KIAA1092 PLCE2

**Immunogen:**

Recombinant Human Inactive phospholipase C-like protein 2 protein (2-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A