

DCTN4 Antibody



PACO42542

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Could have a dual role in dynein targeting and in ACTR1A/Arp1 subunit of dynactin pointed-end capping. Could be involved in ACTR1A pointed-end binding and in additional roles in linking dynein and dynactin to the cortical cytoskeleton.

Gene ID:

DCTN4

Uniprot

Q9UJW0

Synonyms:

Dynactin subunit 4 (Dyn4) (Dynactin subunit p62), DCTN4


Immunogen:

Recombinant Human Dynactin subunit 4 protein (2-460AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: HepG2 whole cell lysate, Jurkat whole cell lysate, HEK293 whole cell lysate, SH-SY5Y whole cell lysate. All lanes: DCTN4 antibody at 3.4 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 53, 46, 54 kDa. Observed band size: 53 kDa.