

HHLA3 Antibody



PACO42574

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

HHLA3

Uniprot

Q9XRX5

Synonyms:

HERV-H LTR-associating protein 3 (Human endogenous retrovirus-H long terminal repeat-associating protein 3), HHLA3

Immunogen:

Recombinant Human HERV-H LTR-associating protein 3 protein (1-121AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A