

IL18 Antibody, Biotin conjugated



PACO42581

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Augments natural killer cell activity in spleen cells and stimulates interferon gamma production in T-helper type I cells By similarity.

Gene ID:

IL18

Uniprot

Q9XSR0

Synonyms:

Interleukin-18 (IL-18) (Interferon gamma-inducing factor) (IFN-gamma-inducing factor) (Interleukin-1 gamma) (IL-1 gamma), IL18, IGIF

Immunogen:

Recombinant Dog Interleukin-18 protein (37-193AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A