

# GINS2 Antibody, HRP conjugated



PACO42587

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

The GINS complex plays an essential role in the initiation of DNA replication, and progression of DNA replication forks. GINS complex seems to bind preferentially to single-stranded DNA.

**Gene ID:**

GINS2

**Uniprot**

Q9Y248

**Synonyms:**

DNA replication complex GINS protein PSF2 (GINS complex subunit 2), GINS2, PSF2

**Immunogen:**

Recombinant Human DNA replication complex GINS protein PSF2 protein (1-185AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A