

# STK38L Antibody, FITC conjugated



PACO42596

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in the regulation of structural processes in differentiating and mature neuronal cells.

**Gene ID:**

STK38L

**Uniprot**

Q9Y2H1

**Synonyms:**

Serine/threonine-protein kinase 38-like (EC 2.7.11.1) (NDR2 protein kinase) (Nuclear Dbf2-related kinase 2), STK38L, KIAA0965 NDR2

**Immunogen:**

Recombinant Human Serine/threonine-protein kinase 38-like protein (212-464AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A