

STK38L Antibody, Biotin conjugated



PACO42597

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the regulation of structural processes in differentiating and mature neuronal cells.

Gene ID:

STK38L

Uniprot

Q9Y2H1

Synonyms:

Serine/threonine-protein kinase 38-like (EC 2.7.11.1) (NDR2 protein kinase) (Nuclear Dbf2-related kinase 2), STK38L, KIAA0965 NDR2

Immunogen:

Recombinant Human Serine/threonine-protein kinase 38-like protein (212-464AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A