

# MRPS28 Antibody, HRP conjugated



PACO42607

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

mitochondrial inner membrane, mitochondrial small ribosomal subunit, mitochondrion, poly(A) RNA binding, mitochondrial translation, mitochondrial translational elongation, mitochondrial translational initiation, mitochondrial translational termination, organelle organization.

**Gene ID:**

MRPS28

**Uniprot**

Q9Y2Q9

**Synonyms:**

28S ribosomal protein S28, mitochondrial (MRP-S28) (S28mt) (28S ribosomal protein S35, mitochondrial) (MRP-S35) (S35mt) (Mitochondrial small ribosomal subunit protein bS1m), MRPS28, MRPS35

**Immunogen:**

Recombinant Human 28S ribosomal protein S28, mitochondrial protein (72-187AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A