

POLDIP2 Antibody



PACO42610

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

mitochondrial nucleoid, mitochondrion.

Gene ID:

POLDIP2

Uniprot

Q9Y2S7

Synonyms:

Polymerase delta-interacting protein 2 (38 kDa DNA polymerase delta interaction protein) (p38), POLDIP2, PDIP38 POLD4


Immunogen:

Recombinant Human Polymerase delta-interacting protein 2 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO42610 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).