

POLDIP2 Antibody, Biotin conjugated



PACO42613

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

mitochondrial nucleoid, mitochondrion.

Gene ID:

POLDIP2

Uniprot

Q9Y2S7

Synonyms:

Polymerase delta-interacting protein 2 (38 kDa DNA polymerase delta interaction protein) (p38), POLDIP2, PDIP38 POLD4

Immunogen:

Recombinant Human Polymerase delta-interacting protein 2 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A