

CARHSP1 Antibody, FITC conjugated



PACO42616

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds mRNA and regulates the stability of target mRNA. Binds single-stranded DNA (in vitro).

Gene ID:

CARHSP1

Uniprot

Q9Y2V2

Synonyms:

Calcium-regulated heat-stable protein 1 (Calcium-regulated heat-stable protein of 24 kDa) (CRHSP-24), CARHSP1

Immunogen:

Recombinant Human Calcium-regulated heat stable protein 1 protein (2-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A